

Product Specifications

ETHANOL (Anhydrous-Alcohol) Biotechnology Grade

Catalog #	<u>Size</u>
IB15720	500ml
IB15721	1L
IB15724	4L

Physical Specifications

CAS# N/A Specific Gravity 0.7964 Max. Formula Weight 46.07 Molecular Formula C₂H₅OH Ethanol 95% Methanol 5% 1% Max. Moisture (KF) Identification (IR) Pass

Molecular Biology Specifications

DNase assay None Detected RNase assay None Detected

Recommended Use

Ethanol is widely used for precipitating nucleic acids. The nucleic acid precipitate, which is formed in the presence of moderate concentrations of monovalent actions, is recovered by centrifugation and redissolved in an appropriate buffer at the desired concentration.

The three major variables associated with this procedure include:

- (1) The temperature at which the precipitate is formed
- (2) Type and concentration of the monovalent cations used in the precipitate mixture.

Most commonly used cations:

- Ammonium acetate (removes the dNTPs therefore do not use as phosorylating a nucleic acid)
- Lithium chloride (used for RNA precipitation)
- Sodium chloride (used for precipitating DNA samples contains SDS)
- Sodium acetate (used for routine DNA and RNA precipitations)
- (3) Time and speed of centrifugation.

Storage

Store at room temperature inside a flame proof cabinet.

Warning

Poison. Flammable liquid. May be harmful or fatal if swallowed. Avoid breathing vapors. Keep away from heat or flame. Do not take internally. See material Safety Data Sheet for additional information.