

Buy Two, Get One FREE


BioDot Serological Pipettes

BioDot Serological Pipettes – Shorter, Ergonomic Design

- 5, 10 & 25ml sizes are 2" shorter for ease of use under the bench top hoods yet long and narrow enough to reach the corners of all flasks
- One piece construction eliminates accidental breakage and minimizes sample hang up
- 2% better accuracy due to construction
- Ascending, descending and negative volume graduations
- Color coded plugged mouth pieces and dispenser boxes for easy identification
- Individually Wrapped, EZ Open, paper/plastic wrap, static free
- Sterile – Certified Nuclease and Nucleic Acid and Pyrogen Free

Catalog #	Size	Length (inches)	Graduation Interval (mL)	Negative Graduation (mL)	Tip I.D. (mm)	Quantity Per/Pkg	Special Promo Price
SPP1-ST-PK	1mL Standard	10.67	0.01	0.3	1.3-1.8	500/Pk	\$73.25
SPP2-ST-PK	2mL Standard	10.67	0.02	0.4	1.8-2.2	500/Pk	\$76.00
SPP5-ST-PK	5mL Standard	11.57	0.10	3.0	1.4-2.0	250/Pk	\$54.50
SPP10-ST-PK	10mL Standard	11.57	0.10	3.0	1.6-2.4	200/Pk	\$46.80
SPP25-ST-PK	25mL Standard	11.67	0.20	2.6	2.3-3.5	100/Pk	\$49.95
SPP50-ST-PK	50mL Standard	14.50	0.50	10.0	3.1-3.3	50/Pk	\$71.40

Serological Promotion valid thru 10-31-11



dot scientific inc.

(800) 878-1785

sales@dotscientific.com

www.dotscientific.com