

# dot scientific inc.

INTRODUCES OUR NEWEST PREMIUM TOUCH® GLOVE

## Premium Touch® Sensitive

Cobalt Blue Nitrile Exam Gloves – Powder Free

“An Excellent Alternative to Latex”

Premium Touch® Sensitive Nitrile Exam Gloves provide superior barrier protection

Feels like Latex – Soft and flexible providing excellent tactile sensitivity

Fits like Latex – Stretchy and comfortable



100% Nitrile – Free of Natural Rubber Proteins  
Thickness 2.5 mil (palm) / 4.0 mil (fingers)  
Textured Finger Tips  
Beaded Cuff - 9" Length  
Class I Medical Device  
Ambidextrous  
**200 Gloves Per Box, 2000 Per Case**

Catalog #	Description	List Price	Introductory Offer
2700-PF	Small	\$160.00	\$120.00
2701-PF	Medium	\$160.00	\$120.00
2702-PF	Large	\$160.00	\$120.00
2703-PF	Extra-Large	\$160.00	\$120.00

For a FREE Sample simply call our customer service representatives at 800-878-1785 or visit us at [www.dotscientific.com](http://www.dotscientific.com)

## dot scientific inc.

(800) 878-1785

[orders@dotscientific.com](mailto:orders@dotscientific.com)

[www.dotscientific.com](http://www.dotscientific.com)

PREMIUM TOUCH® MEDICAL GRADE EXAMINATION GLOVES ARE DESIGNED AND MANUFACTURED TO SPECIFICALLY MEET THE UNIQUE NEEDS OF THE MEDICAL AND RESEARCH ENVIRONMENT. PROVIDING MAXIMUM GRIP, SUPERIOR PROTECTION AND EXCELLENT FIT, THE PREMIUM TOUCH® FAMILY OF GLOVES WILL MEET YOUR RIGOROUS DAILY DEMANDS WITH INCREDIBLE COMFORT.

# Premium Touch®

Blue Nitrile Exam Gloves – Powder Free



100% Nitrile, Free of Natural Rubber Proteins  
 Fully Textured & Odor Free  
 Class I Medical Device  
 5ml Thickness  
 Length – 9 inches  
 Excellent for People with Latex Sensitivities  
 Resistant to a Broad Range of Chemicals  
**100 Gloves Per Box, 1000 Per Case**

Catalog #	Description	List Price	Special Offer
1720-PF	Small	\$105.00	\$72.00
1721-PF	Medium	\$105.00	\$72.00
1722-PF	Large	\$105.00	\$72.00
1723-PF	Extra-Large	\$105.00	\$72.00

# Premium Touch® Gold Grip Plus

Latex Exam Gloves – Powder Free



Ultra-Low Latex Protein Content – 50mg or less  
 6ml Thickness  
 Length – 9.6 inches  
 Class I Medical Device  
 Full Microtexture  
 Excellent Tactile Sensitivity  
**100 Gloves Per Box, 1000 Per Case**



The gloves that labels won't stick to!

Catalog #	Description	List Price	Special Offer
1350-XSPF	Extra Small	\$98.00	\$72.00
1350-PF	Small	\$98.00	\$72.00
1351-PF	Medium	\$98.00	\$72.00
1352-PF	Large	\$98.00	\$72.00
1353-PF	Extra-Large	\$98.00	\$72.00