

MVE CryoShipper CT-50

The MVE CryoShipper CT-50 is specifically engineered to ship one to two 25 mL blood bags. Our unique rack design that holds one to two blood bags helps limit sample movement during transit.

Features include:

- Cassette Retrieval
- Secure shipping for one to two 25 mL blood bags
- Datalogger option is available
- Charges in 2 hours
- Safe for international shipping; IATA Dangerous Goods Regulations Exempt (Special Provision A152)



MVE CryoShipper CT-50

Maximum Storage Capacity

Number of Blood Bags 25 mL	2
----------------------------	---

Performance

LN2 Capacity (min) L est.	5
---------------------------	---

Static Evaporation Rate* L/day	0.5
--------------------------------	-----

Static Holding Time* days	10
---------------------------	----

Unit Dimensions

Neck Opening in (mm)	3.82 (97)
----------------------	-----------

Overall Height in (mm)	15 (380)
------------------------	----------

Outside Diameter in (mm)	14.5 (368)
--------------------------	------------

Weight Empty lb. (kg)	17.6 (8)
-----------------------	----------

Weight Charged vapor lb. (kg)	26.45 (12)
-------------------------------	------------

* Static evaporation rate and static holding time are nominal. Actual rate and holding time will be affected by the nature of container use, atmospheric conditions, and manufacturing tolerances.

CT = Cell Therapy

TWO Year Tank Warranty

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU.

Visit www.chartbiomed.com for more information.