

MVE CryShipper CT-250

The MVE CryShipper CT-250 is engineered specifically to be able to ship one to two 250 mL blood bags. Our unique rack design helps limit sample movement during transit. This cryoshipper holds one cassette retrieval rack designed to hold one or two 250 mL cassettes.

Features include:

- CT250 Frame puller
- Secure shipping for one to two 250 mL blood bags
- Charges in under 2 hours
- Safe for international shipping; IATA Dangerous Goods Regulations Exempt (Special Provision A152)



MVE CryShipper CT-250

Maximum Storage Capacity

Number of Blood Bags 250 mL	2
-----------------------------	---

Performance

LN2 Capacity (min) L	8
----------------------	---

Static Evaporation Rate* L/day	0.8
--------------------------------	-----

Static Hold Time* days	10
------------------------	----

Unit Dimensions

Neck Opening in (mm)	6.5 (165)
----------------------	-----------

Overall Height in (mm)	19.8 (502)
------------------------	------------

Outside Diameter in (mm)	16.0 (406)
--------------------------	------------

Weight Empty lb. (kg)	24.3 (11)
-----------------------	-----------

Weight Charged vapor lb. (kg)	38.6 (17.5)
-------------------------------	-------------

Frame Puller Height in (mm)	15.2 (385)
-----------------------------	------------

Frame Puller Width in (mm)	6.3 (161)
----------------------------	-----------

* Static evaporation rate and static holding time are nominal. Actual rate and holding time will be affected by the nature of container use, atmospheric conditions, and manufacturing tolerances.

CT = Cell Therapy

TWO Year Tank Warranty

Conforms to MDD 93/42/EEC, the Medical Device Directive for the EU.

Visit www.chartbiomed.com for more information.

Copyright © 2016 Chart Industries Ref ML-A10032 B

Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.