

# ProPette™ Controller

## ELECTRONIC PIPETTE CONTROLLER

### WITH QUICKSTAND™



- **With QuickStand Retractable Stand**
- **Variable Speed Control**
- **Cordless and rechargeable**
- **8 hours of continuous operation**
- **For pipettes from 1ml to 100ml**

ProPette™ Controllers from MTC Bio™ are lightweight (only 6.3oz/170g) and designed to fit comfortably in either the right or the left hand. The natural feel of the ergonomic body design, the computer balanced weight distribution and the low pressure fingertip control buttons...all combine to provide a more comfortable pipetting experience.

A long-life rechargeable lithium (environmentally friendly) battery keeps the controller running for up to 8 hours or 2500 pipettings on a single two hour charge. ProPette's easy-fit silicone nosepiece accepts pipette sizes from 1ml to 100ml. With its light weight advanced ergonomics, continuously variable electronic speed adjustment and extended lithium battery life, the ProPette Pipette Controller is the logical choice for comfort, convenience and value.

**MTC Bio**

2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 [www.mtcbio.com](http://www.mtcbio.com)

Effective March 1, 2013. Subject to change without notice.

### Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Price</u>
P6080	ProPette XL Pipette Controller; includes charger and US type plug (100-240V for worldwide use) benchtop stand, extra filter and replaceable lithium battery (5V/700mA)	\$279.00
P6080-I	ProPette XL Pipette Controller; as above, but with Euro-style 2 prong plug	\$279.00
P6080-BA	Replacement battery, lithium, 5.0V, 700mA	\$29.00
P6080-FI	Replacement filter, 0.45um	\$5.00