

2 Liter General Purpose Water Bath



WBE02

Key Specifications

Description	Gen. Purpose Water Bath, Digital, 2 L
Display	Full-color TFT
Languages Supported	5: French, German, Spanish, Chinese, English
Working Temperature Range °C	Ambient +5° to 100°
Temperature Stability °C	±0.1°
Temperature Calibration Capability	1-point
Working Access (L x W x D) (cm)	9.9 x 10.9 x 15.2
Flammability Class (DIN 12876-1)	I (NFL)

Key Features

Common Specifications

Working Temperature Range °F	Ambient +10° to 212°
Working Temperature Range °C	Ambient +5° to 100°
Reservoir Capacity (gallons)	0.53
Reservoir Capacity (liters)	2
Reservoir/Tank Material	Stainless Steel
Reservoir Cover	Hinged Gable
Working Access (L x W x D) (inches)	3.9 x 4.3 x 6

Working Access (L x W x D) (cm)	9.9 x 10.9 x 15.2
Temperature Stability °F	±0.2°
Temperature Stability °C	±0.1°
Display	Full-color TFT
Flammability Class (DIN 12876-1)	I (NFL)
Over-Temperature Protection / Failsafe Heater Control	Yes
Maximum Ambient Temperature °F	104°
Maximum Ambient Temperature °C	40°
Overall Dimensions (L x W x H) (inches)	8.8 x 8.5 x 13.1
Overall Dimensions (L x W x H) (cm)	22.2 x 21.6 x 33.3
Shipping Weight (pounds)	9.0
Shipping Weight (kilograms)	4.1
Catalog Page Number	86-87

50 Hz Only

Part Number	WBE02A12E
Included Hardware	Sample Tray
Heater Wattage	120
Electrical Requirements (VAC/Hz/Ph/A)	240/50/1/0.5
Regulatory Approvals	CE



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 05/20/2020 - 11:36): <https://www.polyscience.com/general-purpose-water-baths/2-liter-general-purpose-water-bath>