

## 10 Liter General Purpose Water Bath



WBE10

## Key Specifications

<b>Description</b>	Gen. Purpose Water Bath, Digital, 10 L
<b>Display</b>	Full-color TFT
<b>Languages Supported</b>	5: French, German, Spanish, Chinese, English
<b>Working Temperature Range °C</b>	Ambient +5° to 100°
<b>Temperature Stability °C</b>	±0.1°
<b>Temperature Calibration Capability</b>	1-point
<b>Working Access (L x W x D) (cm)</b>	26.9 x 29.5 x 15.2
<b>Flammability Class (DIN 12876-1)</b>	I (NFL)

## Key Features

## Common Specifications

<b>Working Temperature Range °F</b>	Ambient +10° to 212°
<b>Working Temperature Range °C</b>	Ambient +5° to 100°
<b>Reservoir Capacity (gallons)</b>	2.64
<b>Reservoir Capacity (liters)</b>	10
<b>Reservoir/Tank Material</b>	Stainless Steel
<b>Reservoir Cover</b>	Hinged Gable
<b>Working Access (L x W x D) (inches)</b>	10.6 x 11.6 x 6

<b>Working Access (L x W x D) (cm)</b>	26.9 x 29.5 x 15.2
<b>Temperature Stability °F</b>	±0.2°
<b>Temperature Stability °C</b>	±0.1°
<b>Display</b>	Full-color TFT
<b>Flammability Class (DIN 12876-1)</b>	I (NFL)
<b>Over-Temperature Protection / Failsafe Heater Control</b>	Yes
<b>Reservoir Drain</b>	Yes
<b>Maximum Ambient Temperature °F</b>	104°
<b>Maximum Ambient Temperature °C</b>	40°
<b>Overall Dimensions (L x W x H) (inches)</b>	16.3 x 14.4 x 15.0
<b>Overall Dimensions (L x W x H) (cm)</b>	41.3 x 36.5 x 38.1
<b>Shipping Weight (pounds)</b>	22.0
<b>Shipping Weight (kilograms)</b>	10.0
<b>Catalog Page Number</b>	86-87

## 50 Hz Only

<b>Part Number</b>	WBE10A12E
<b>Included Hardware</b>	Sample Tray
<b>Heater Wattage</b>	1000
<b>Electrical Requirements (VAC/Hz/Ph/A)</b>	240/50/1/4.5
<b>Regulatory Approvals</b>	CE



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

**Phone:** +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

\*/

**Source URL (modified on 05/20/2020 - 11:42):** <https://www.polyscience.com/general-purpose-water-baths/10-liter-general-purpose-water-bath>