

20 Liter General Purpose Water Bath



WBE20

Key Specifications

| | |
|---|--|
| Description | Gen. Purpose Water Bath, Digital, 20 L |
| Display | Full-color TFT |
| Languages Supported | 5: French, German, Spanish, Chinese, English |
| Working Temperature Range °C | Ambient +5° to 100° |
| Temperature Stability °C | ±0.1° |
| Temperature Calibration Capability | 1-point |
| Working Access (L x W x D) (cm) | 24.1 x 43.2 x 15.2 |
| Flammability Class (DIN 12876-1) | I (NFL) |

Key Features

Common Specifications

| | |
|--|----------------------|
| Working Temperature Range °F | Ambient +10° to 212° |
| Working Temperature Range °C | Ambient +5° to 100° |
| Reservoir Capacity (gallons) | 5.28 |
| Reservoir Capacity (liters) | 20 |
| Reservoir/Tank Material | Stainless Steel |
| Reservoir Cover | Hinged Gable |
| Working Access (L x W x D) (inches) | 9.5 x 17 x 6 |

| | |
|--|--------------------|
| Working Access (L x W x D) (cm) | 24.1 x 43.2 x 15.2 |
| Temperature Stability °F | ±0.2° |
| Temperature Stability °C | ±0.1° |
| Display | Full-color TFT |
| Flammability Class (DIN 12876-1) | I (NFL) |
| Over-Temperature Protection / Failsafe Heater Control | Yes |
| Reservoir Drain | Yes |
| Maximum Ambient Temperature °F | 104° |
| Maximum Ambient Temperature °C | 40° |
| Overall Dimensions (L x W x H) (inches) | 15.1 x 22.5 x 15.1 |
| Overall Dimensions (L x W x H) (cm) | 38.4 x 57.2 x 38.4 |
| Shipping Weight (pounds) | 25.0 |
| Shipping Weight (kilograms) | 11.3 |
| Catalog Page Number | 86-87 |

50 Hz Only

| | |
|--|-------------|
| Part Number | WBE20A12E |
| Included Hardware | Sample Tray |
| Heater Wattage | 1400 |
| Electrical Requirements (VAC/Hz/Ph/A) | 240/50/1/6 |
| Regulatory Approvals | CE |



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 05/21/2020 - 15:30): <https://www.polyscience.com/general-purpose-water-baths/20-liter-general-purpose-water-bath>