

28 Liter General Purpose Water Bath



WBE28

Key Specifications

Description	Gen. Purpose Water Bath, Digital, 28 L
Display	Full-color TFT
Languages Supported	5: French, German, Spanish, Chinese, English
Working Temperature Range °C	Ambient +5° to 100°
Temperature Stability °C	±0.1°
Temperature Calibration Capability	1-point
Working Access (L x W x D) (cm)	24.1 x 43.2 x 20.3
Flammability Class (DIN 12876-1)	I (NFL)

Key Features

Common Specifications

Working Temperature Range °F	Ambient +10° to 212°
Working Temperature Range °C	Ambient +5° to 100°
Reservoir Capacity (gallons)	7.4
Reservoir Capacity (liters)	28
Reservoir/Tank Material	Stainless Steel
Reservoir Cover	Hinged Gable
Working Access (L x W x D) (inches)	9.5 x 17 x 8

Working Access (L x W x D) (cm)	24.1 x 43.2 x 20.3
Temperature Stability °F	±0.2°
Temperature Stability °C	±0.1°
Display	Full-color TFT
Flammability Class (DIN 12876-1)	I (NFL)
Over-Temperature Protection / Failsafe Heater Control	Yes
Reservoir Drain	Yes
Maximum Ambient Temperature °F	104°
Maximum Ambient Temperature °C	40°
Overall Dimensions (L x W x H) (inches)	15.1 x 22.5 x 17.0
Overall Dimensions (L x W x H) (cm)	38.4 x 57.2 x 43.2
Shipping Weight (pounds)	28.0
Shipping Weight (kilograms)	12.7
Catalog Page Number	86-87

50 Hz Only

Part Number	WBE28A12E
Included Hardware	Sample Tray
Heater Wattage	1400
Electrical Requirements (VAC/Hz/Ph/A)	240/50/1/6
Regulatory Approvals	CE



© 2015-2020 **PolyScience** 6600 W. Touhy Avenue, Niles, Illinois 60714-4516 USA

Phone: +1(847) 647-0611 or +1(800) 229-7569 **Fax:** +1(847) 647-1155

*/

Source URL (modified on 05/20/2020 - 11:44): <https://www.polyscience.com/general-purpose-water-baths/28-liter-general-purpose-water-bath>