Product Name:	ABT-HC-SLS-16
Product Name: Description:	ABT-HC-SLS-16 Standard Solid Door Laboratory Refrigerator 16 Cu. Ft. capacity One swing solid door, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: $1^{\circ}C - 10^{\circ}C$ Three adjustable shelves (adjustable in ½" increments) Forced draft air circulation Heavy duty door handle Keyed door lock Magnetic door gasket for positive seal LED interior lights are safety shielded and switch controlled Swivel Casters ( $3^{1/2}$ ") Probe access port ( $3/8$ ") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 25" W x 29 <sup>3</sup> / <sub>4</sub> " D x 79" H with casters; 75 <sup>1</sup> / <sub>2</sub> " without casters Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 277 lbs.