Product Name:	ABT-HC-SLS-23
Description:	Standard Solid Door Laboratory Refrigerator
	23 Cu. Ft. capacity
	One swing solid door, self-closing, right hinged, not reversible
	Microprocessor temperature controller
	Digital temperature display
	Audible and visual high and low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Remote alarm contacts
	Adjustable operating temperature range: 1°C – 10°C
	Keyed door lock
	Magnetic door gasket for positive seal
	Heavy duty door handle
	Four adjustable shelves (adjustable in ½" increments) Forced draft air circulation
	LED interior lights are safety shielded and switch controlled
	Swivel Casters (3 ½")
	Probe access port (3/8")
	White powder-coated steel exterior and interior
	Hydrocarbon, natural refrigerant (R290)
	One year parts and labor warranty, plus an additional four year
	compressor parts warranty
	Exterior dimensions: 26 1/8" W x 34 3/4" D x 81 3/4" H with casters; 78
	1/4" H without casters
	Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP
	ETL listed
	Shipping weight: 250 lbs.