Product Name:	ABT-HC-CS-33
Description:	Standard Glass Door Chromatography Refrigerator 33 Cu. Ft. Capacity Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)
	Microprocessor temperature controller
	Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts
	Adjustable operating temperature: 1°C to 10°C Eight adjustable shelves (adjustable in 1/2" increments) Forced draft air circulation
	Keyed door lock Magnetic door gasket for positive seal
	Integrated door handles Swivel casters (3 ½") LED interior lights are safety shielded and switch controlled
	One GFCI protected duplex electrical outlet One 2" access port with spring loaded self-closing, gasketed metal
	cover Probe access port (3/8")
	Powder coated steel white interior and exterior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290)
	One year parts and labor warranty, plus an additional four year compressor parts warranty
	Exterior dimensions: 40 3/8" W* x 33 1/4" D x 81 3/4" H with casters; 78 1/4" H without casters. *Width w/access port closed; 44 3/8" width with access port open
	Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/2 HP
	ETL listed Shipping weight: 419 lbs.