

Product Name:	ABT-HC-CS-47
Description:	<p>Standard Glass Door Chromatography Refrigerator  47 Cu. Ft. Capacity  Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature: 1°C to 10°C  Eight adjustable shelves (adjustable in ½" increments)  Forced draft air circulation  Keyed door lock  Magnetic door gasket for positive seal  Integrated door handles  Swivel casters (3 ½")  LED interior lights are safety shielded and switch controlled  One GFCI protected duplex electrical outlet  One 2" access port with spring loaded self-closing, gasketed metal cover  Probe access port (3/8")  Powder coated steel white interior and exterior  High density urethane foam cabinet insulation  Hydrocarbon, natural refrigerant (R290)  One year parts and labor warranty, plus an additional four year compressor parts warranty  Exterior dimensions: 54 7/8" W* x 33 ¼" D x 81 ¾" H with casters; 78 ¼" H without casters. *Width w/access port closed; 58 ¾"W width with access port open  Cycle defrost  115V, 60 Hz, 4.5 Amps, 1/2 HP  ETL listed  Shipping weight: 508 lbs.  Note: This unit requires 2 outlets</p>