Product Name:	PH-ART-HC-RFC9
T TOGGET NAME.	
Product Name: Description:	Premier Pharmacy/Vaccine Refrigerator & Freezer Combination 9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer) Two exterior solid swing doors right hinged* Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C25°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm Audible and visual high and low temperature alarms, refrigerator and freezer Remote alarm contacts, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall
	2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being
	Keyed door locks Magnetic door gasket for positive seal Recessed door handle, freezer; heavy duty door handle, refrigerator Leveling legs Stacking kit LED interior lights Two probe access ports (3/8"); one, refrigerator/one, freezer
	White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 23 ¾"W x 24 ½"D x 66 ¾"H Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer UL/C-UL listed
	Shipping weight: 269 lbs. *Please note this unit requires 2 outlets. *This product is available left hinged, product number PH-ABT-HC-RFC9-LH