

**New!**

# ClearSeal™

## Screw-Cap Microtubes



- Ultra-clear polypropylene
- Packaged in *resealable* bags of 50
- Leak-proof O-ring seal
- Sterile, RNase, DNase & Pyrogen free
- White, printed graduations

ClearSeal™ screw-cap microtubes come sterile, in a variety of sizes, cap configurations, graduation options, and even a loop-cap option for the 2.0mL version. Made from ultra-clear polypropylene, these screw-cap microtubes include assembled caps with silicone O-rings that assure a leakproof seal from -90°C to +121°C. All ClearSeal™ screw cap microtubes are sterile, RNase, DNase, & pyrogen free and may be centrifuged at up to 22,000xg\*.

\*Refers to liquids with a specific gravity of 1.0, using appropriate rotors, cushions, and adapters.



Packaged in *resealable* bags of 50

Cat. No.	Description	Quantity
C3205-C	ClearSeal™ 0.5mL (10.2 x 36.5mm), sterile, w/ O-ring & attached caps, non-graduated, conical bottom, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3205-S	ClearSeal™ 0.5mL (13 x 47.3mm), sterile, w/ O-ring & attached caps, non-graduated, self-standing, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3215-CG	ClearSeal™ 1.5mL (12.8 x 45.3mm), sterile, w/ O-ring & attached caps, non-graduated, conical bottom, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3215-SG	ClearSeal™ 1.5mL (12.8 x 47.6mm), sterile, w/ O-ring & attached caps, printed graduations, self-standing, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3220-CG	ClearSeal™ 2.0mL (12.9 x 46.5mm), sterile, w/ O-ring & attached caps, printed graduations, conical bottom, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3220-SG	ClearSeal™ 2.0mL (13.1 x 46.8mm), sterile, w/ O-ring & attached caps, printed graduations, self-standing, 50 per sterile resealable bag, 20 bags per case	1000/cs
C3220-SGL	ClearSeal™ 2.0mL (13 x 49.2mm), sterile, w/ O-ring & attached loop-caps, printed graduations, self-standing, 50 per sterile resealable bag, 20 bags per case	1000/cs

