

Premium LT Pipette

NEW

70% less pipetting forces



Glass-made counter window

More durable than plastic and easy to read.



Lock mechanism

Stronger material used in 4 notches improves lock strength, which minimizes the lever drifting and wearing.



Plunger

Ceramic plunger is used to prevent corrosion by solvents. (100µL & up)



Body

Material used has high tolerance to solvents and to drop impact.



Nozzle

Nozzle top is highly abrasion tolerant for 2µL through 1000µL.



Hyper blower system

Hyper blower system improves dispensing efficiency, especially with 2µL and 10µL.

- Light Weight & Ergonomic Shape
- Noncorrosive Ceramic Plunger 100ul and up
- Fully Autoclavable & UV Resistant
- Volume Lock
- Easy Calibration
- Hyper Blower System
- Universal Nozzle
Specialty tips are not required
- Also Available in White
Add -W to part# or DS-NLTW-4PK

DOT Scientific Premium LT Pipettes 5 Year Warranty

Catalog#	Description	Size	List	Special Offer
DS-2-NLT	Single Channel Pink	0.1 - 2ul Pipette	\$ 450.00	\$ 230.00
DS-10-NLT		0.5 - 10ul Pipette	\$ 434.00	\$ 215.00
DS-20-NLT		2 - 20ul Pipette	\$ 434.00	\$ 215.00
DS-100-NLT		10 - 100ul Pipette	\$ 434.00	\$ 215.00
DS-200-NLT		20 - 200ul Pipette	\$ 434.00	\$ 215.00
DS-1000-NLT		100 - 1000ul Pipette	\$ 434.00	\$ 215.00



+



= Savings 4 You!

\$ 790.00

- Includes #DS-MLT-STD2 6 Place Rotary Stand
- Custom 4 Packs Available for \$800 - choose from all 6 models

DOT Scientific Premium LT Pipettes 4 Pack and Stand

Catalog#	Description	List	Special Offer
DS-NLT-4PK	Includes 10, 20, 200 & 1000ul Vols. & Stand	\$ 980.00	\$ 790.00

DOT Scientific, Inc. • 4165 Lippincott • Burton, MI 48519 • Toll Free 800.878.1785

Fax 810.744.2865 • orders@dotscientific.com • www.dotscientific.com

Expires 12-31-2022

Premium EXII Pipette



- Fully Autoclavable
Ideal for tissue culture & RNA
- UV Resistant
- Volume Lock
Set volume will not drift
(Good for beginners)
- Noncorrosive Ceramic Plunger
200ul & up
- Replaceable Parts
Long service life
- Thermally Insulated Handle
Minimize changes in accuracy due to hand temperature
- Universal Nozzle
Specialty tips are not required
Nozzle filtered on 1000ul & up
- Ask About Premium EXPlusII
Solvent, chemical resistant models

DOT Scientific Premium EXII Pipettes 2 Year Warranty

Catalog#	Description	Size	List	Special Offer
DS-2-NPX2	Single Channel Blue	0.1 - 2ul Pipette	\$ 346.97	\$ 185.00
DS-10-NPX2		0.5 - 10ul Pipette	\$ 336.97	\$ 175.00
DS-20-NPX2		2 - 20ul Pipette	\$ 336.97	\$ 175.00
DS-100-NPX2		10 - 100ul Pipette	\$ 336.97	\$ 175.00
DS-200-NPX2		20 - 200ul Pipette	\$ 336.97	\$ 175.00
DS-1000-NPX2		100 - 1000ul Pipette	\$ 336.97	\$ 175.00
DS-5000-NPX2		1000 - 5000ul Pipette	\$356.97	\$ 195.00
DS-10000-NPX2		1000 - 10000ul Pipette	\$356.97	\$ 195.00



+ = Savings 4 You!
\$ 665.00

- Includes #DS-MLT-STD2
6 Place Rotary Stand
- Custom 4 Packs Available
for \$680 - choose from all
8 models

DOT Scientific EXII Pipettes 4 Pack and Stand

Catalog#	Description	List	Special Offer
DS-NPX2-4PK	Includes 10, 20, 200 & 1000ul Vols. & Stand	\$ 980.00	\$ 665.00

DOT Scientific, Inc. • 4165 Lippincott • Burton, MI 48519 • Toll Free 800.878.1785
Fax 810.744.2865 • orders@dotscientific.com • www.dotscientific.com

Expires 12-31-2022