

EBA 200 | S

BLOOD / URINE / PEDIATRIC / PPP / COAG.

Compact performance with built-in versatility for practical centrifugation.

The EBA 200 is a compact clinical centrifuge offering performance and safety for small sample volumes. The EBA 200 features a built-in, versatile, 8-place, fixed-angle rotor and adjustable LCD display for convenient and quick operation. Available in low and high speed (up to 8,000 RPM) models.

High-performance for quick turnaround of time sensitive samples.

The higher performance EBA 200 S version delivers an RCF of 6,153 making it ideal for STAT chemistry, PPP and Coag. Achieve consistent platelet counts in just minutes for Platelet Poor Plasma (PPP).

Adjustable parameters for flexible application control.

The EBA 200 and EBA 200 S feature a flexible control panel allowing you to adjust your application parameters. Adjust time, speed (RPM), RCF and braking with the easy to use LCD display.

Visual and audible end-of-cycle notification for a productive workflow.

Customizable and independent control of visual and/or audio alarms can be set for end of cycle notification. This provides an across-the-room indication that your unit is ready for the next operation.



MODEL	EBA 200	EBA 200 S
Max. Capacity	8 x 15 ml	8 x 15 ml
Max. RPM / RCF	6,000 / 3,461	8,000 / 6,153
Max. Blood Tubes*	8 (up to 10 ml)	8 (up to 10 ml)
Max. Conical Tubes*	4 (15 ml)	4 (15 ml)
Temperature Control	Ambient	Ambient
Dims. (H x W x D)	9 x 10 x 14 in.	9 x 10 x 14 in.
Warranty†	5 yr.	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more at: hettweb.com/eba200