Product Name:	ABT-MFP-30-TS
Description:	TempLog Premier Manual Defrost Laboratory Freezer
	30 Cu. Ft. capacity
	One solid swing door, right hinged, not reversible
	Five inner "fast-freeze compartments" with doors for temperature maintenance
	and recovery
	Microprocessor Temperature Controller
	Adjustable Operating Temperature Range: 15°-125°C
	Touch Screen Digital Temperature Display and Alarm Module
	Password-protected Control Panel Sample simulated digital temperature display; .1° Resolution, adjustable
	temperature display in Celsius or Fahrenheit
	Active Temperature Scale Bar Display
	Resettable 1-point Min/Max History
	Digital 24-hour Temperature Chart
	Temperature Data Logging
	USB Port for Data Transfer and Software Updates
	Real-time Operational Display Clock
	4-20mA Output
	Electronic Troubleshooting Guide
	Electronic QR Code to Access Owner's Manual & Technical Support Contact
	Information
	Adjustable Alarm Volume
	Audible and Visual High/Low Temperature Alarms
	Sensor Error Alarm
	Power Failure Alarm
	Door Ajar Alarm Low Battery Alarm
	Operation Indicator
	Alarm Validation (audible, visual, and remote alarm)
	Alarm Mute with Ring Back
	Remote Alarm Contacts
	Four fixed shelves
	Keyed door lock
	Magnetic door gasket for positive seal
	Stainless steel door handle
	Caster platform (adds 4" to height)
	One probe access port (3/8")
	White textured cabinet and door
	High density urethane foam cabinet and door insulation
	CFC-free refrigerant (R404a), foam insulation and packaging
	Three year parts and labor warranty, plus additional two year compressor parts
	warranty Exterior dimensions: 35 3/"W x 3.4 1/"D x 76 1/"H with caster platform: 72 1/"H
	Exterior dimensions: 35 ¾"W x 34 ¼"D x 76 ½"H with caster platform; 72 ½"H without caster platform
	Manual defrost
	115V, 60 Hz, 6 Amps, 1/3 HP
	UL/C-UL listed
	Shipping weight: 462 lbs.
L	- U - U