

# SAVE!!!

**The LabMart®**  
Serving Science Since 1850

# Transfer Pipets



- Low density PE; non-toxic and inert
- Attached bulb to avoid cross-contamination
- One-step biohazard disposal
- Can be frozen in liquid nitrogen
- May be used in procedures with beakers, dishes, flasks, slides, vials, test tubes or roller bottles
- Ideal for use with volumetric flasks, multiple well plates, centrifuge tubes or cryogenic vials

Purchase 2 Cases and  
Receive a FREE Beaker  
Mug



		<b>Box WAS</b>	<b>Box NOW</b>	<b>Case WAS</b>	<b>Case NOW</b>
<b>M373570</b>	Transfer Pipets Standard 7.7ml without Graduations	\$85.37	<b>\$21.18/500</b>	\$768.30	<b>\$241.96/5000</b>
<b>M373570-S</b>	Sterile 7.0ml Transfer Pipet 3.5ml Draw	174.47	<b>97.83/400</b>	875.00	<b>559.00/4000</b>
<b>M373571</b>	Transfer Pipet Blood Bank 5ml, 6 in.	66.92	<b>18.76/500</b>	566.80	<b>214.28/5000</b>
<b>M373572</b>	Transfer Pipet Graduated 1ml Large Bulb 5.8ml Cap.	74.01	<b>37.69/500</b>	626.38	<b>215.34/5000</b>
<b>M373572-S</b>	Sterile Transfer Pipet Graduated 1ml Total Cap. 5.8ml	207.44	<b>97.30/400</b>	865.00	<b>556.00/4000</b>
<b>M373572-2S</b>	Sterile Graduated Transfer Pipet	146.00	<b>37.49/400</b>	1,167.97	<b>214.22/4000</b>
<b>M373573</b>	3ml Graduated Transfer Pipet	66.92	<b>34.51/500</b>	1,167.97	<b>197.16/5000</b>
<b>M373574</b>	2ml Graduated Transfer Pipet	82.09	<b>29.44/500</b>	566.80	<b>168.22/5000</b>
<b>M373575</b>	Fine Tip Transfer Pipet 5.8ml Cap.	129.51	<b>82.22/500</b>	634.66	<b>469.78/5000</b>
<b>M373578</b>	Transfer Pipet 4.6ml Cap. 3.3ml Draw Narrow Stem Standard Tip	66.10	<b>40.43/500</b>	1,097.07	<b>231.00/5000</b>
<b>M373579</b>	Transfer Pipet 6ml Cap. 2.3ml Draw Extra Long 9 in.	123.90	<b>59.02/500</b>	1,115.05	<b>337.22/5000</b>
<b>M373579-S</b>	Sterile Extra Long Transfer Pipet	225.66	<b>113.40/400</b>	980.00	<b>648.00/4000</b>
<b>M773570</b>	Transfer Pipet 23ml Cap. 7.3ml Draw Extra Long 12 in.	110.09	<b>36.00/100</b>	990.79	<b>205.66/1000</b>

**HURRY Offer ends 3/31/24**