



ReadyStrain[™] Cell Strainer Kits

- Includes Strainer, 50mL Tube, and Tube Cap
- Preassembled, sterile and ready to use
- Packaged in easy-open “tear-strip” bags
- Ready to open directly under the hood
- Gamma sterilized to SAL 10⁻⁶



Watch the video demonstration!



Each sterile kit contains a SureStrain[™] strainer, pre-attached to a 50mL tube, with a separate screw-cap ready for assembly, processing and storage after the filtration is complete.

Kits are supplied in easy open “tear-strip” bags and gamma sterilized.



ReadyStrain[™] is the only product of its kind, providing all necessary parts for removing clumps and debris to isolate single cell suspensions. ReadyStrain[™] supplies all essential components, combining them into a single, preassembled, ready to use, sterile packet. These kits eliminate the contamination risks associated with opening multiple packages and assembling separate pieces (tubes, strainers and caps). The individual sterile packets mean researchers are no longer forced to expose entire bags of 50mL tubes and caps in order to perform one filtration. Instead, they can open a single packet directly under the hood.

CAT. NO.	DESCRIPTION	QTY	PRICE
C5040	Cell Strainer Kit: Blue 40µm Strainer, 50mL Tube and Screw-Cap, Sterile	50	\$166.80
C5070	Cell Strainer Kit: White 70µm Strainer, 50mL Tube and Screw-Cap, Sterile	50	\$166.80
C5100	Cell Strainer Kit: Yellow 100µm Strainer, 50mL Tube and Screw-Cap, Sterile	50	\$166.80