

SAVE BIG WITH LABNET SPECIAL OFFERS

# LabValues 2024

From July 1 to December 31, 2024

## Labnet AccuSeal™ Semi-automated Plate Sealer

Offering complete versatility, the AccuSeal will accept a full range of plates for PCR, assay or storage applications, and can be used for standard height to deep well plates.

Sealing parameters are set and displayed via the user-friendly control panel, and the sealing operation is automated to generate consistent results.

- ▶ Adjustable sealing temperature: 100°C to 190°C
- ▶ Precise temperature, timing, and pressure for consistent sealing

Cat. No.	Description	List Price	Special Price
PS1000	AccuSeal semi-automated plate sealer, includes adapters for standard and deep well microplates, 120V	\$4,800.00	<b>\$2,400.00</b>



## Labnet MultiGene™ OptiMax

The MultiGene OptiMax Thermal Cycler delivers advanced speed and features while providing premium performance at an affordable price.

The OptiMax comes with standard built in 2-step, temperature optimization, touchdown and time increment protocols that are easy to adjust to meet your cycling needs. You can select lid temperature or turn it off depending on your needs. The faster ramp rates allow more work to be done in a work day.

- ▶ Unique flexible programming with FlexTemp™ technology
- ▶ Fast ramp rates 5°C heating / 3.5°C cooling
- ▶ Simple user interface

Cat. No.	Description	List Price	Special Price
TC9610	MultiGene OptiMax with 96-well block, 120V	\$4,118.00	<b>\$2,882.60</b>

## Labnet Enduro™ Power Supply

The Enduro Power Supplies are ideal for DNA, RNA and protein electrophoresis. They provide a high degree of sophistication and versatility yet are extremely easy to use and simple to program.

- ▶ Dependable performance, intelligent control
- ▶ Compact, lightweight and stackable

Cat. No.	Description	List Price	Special Price
E0303	Enduro 300V power supply	\$430.00	<b>\$301.00</b>
E0203	Enduro 250V power supply	\$1,184.14	<b>\$828.90</b>

ONLY  
**\$301**  
(300V)



2LV-3.24

## Buy Labnet AccuBlock™ Digital Dry Bath and Get Block for FREE



2LV-4.24



GET  
1 OR 2  
BLOCKS  
FREE

Purchase a Labnet AccuBlock Digital Dry Bath (single block) and get one (1) block at no additional cost.

Purchase a Labnet AccuBlock Digital Dry Bath (dual block) and get two (2) blocks at no additional cost.

Cat. No.	Description
D1301	Digital dry bath, single block, 120V
D1302	Digital dry bath, dual block, 120V

### Get One or Two Blocks for FREE

Cat. No.	Description
D1101	Solid block for machining (no holes)
D1102	Single block, 48 x 0.2 mL PCR tubes or 6 x 0.2 mL
D1102A	Single block, 24 x 2.0 mL tubes
D1105	Single block, 24 x 0.5 mL tubes
D1105A	Single block, 24 x 1.5 mL tubes
D1106	Single block, 35 x 6 mm tubes
D1110	Single block, 20 x 10 mm tubes
D1112	Single block, 20 x 12 mm tubes

Cat. No.	Description
D1113	Single block, 20 x 13 mm tubes
D1115-TALL	Single block, 12 x 15 mL centrifuge tubes
D1116	Single block, 12 x 15 or 16 mm tubes
D1117	Single block, 12 x 17 mm tubes
D1120	Single block, 6 x 20 mm tubes
D1125	Single block, 6 x 25 mm tubes
D1150-TALL	Single block, 5 x 50 mL conical tubes
D1196-PCR	Single block, 96-well PCR plate, skirted or non-skirted
D1296	Dual block, 96-well microplate or 4 slides
D1296-PCR	Dual block, 96-well PCR plate, skirted or non-skirted
D12384	Dual block, 384-well PCR plate (for dual block unit only)


## Labnet Enduro™ GDS II Gel Documentation Systems

Enduro GDS II and GDS Touch II are used for visualization of DNA and protein stained within a gel and provide photo documentation of these gels in a laboratory setting.

The GDS II comes in either a 302 nm or 365 nm version and runs on a Windows-based PC (not included). The GDS Touch II comes in either a 302 nm or 365 nm wavelength version with 470 nm epi-blue lights for use with “safe dyes” that eliminate the need for ethidium bromide.

- Region of interest function allows optimal exposure in less than 30 seconds. No more multiple exposures to get good images.

Cat. No.	Description	List Price	Special Price
GDS2-1302	Enduro GDS II with 302 nm UV	\$6,073.92	<b>\$4,859.14</b>
GDS2-1365	Enduro GDS II with 365 nm UV	\$6,073.92	<b>\$4,859.14</b>
GDST2-1302	Enduro GDS Touch II with epi-blue and 302 nm UV	\$9,874.08	<b>\$7,899.26</b>



2LV-5.24



2LV-6.24

## Labnet Enduro™ UV Transilluminators

The unique hinge design and frame on the UV shield prevents the lid from cracking and hinge damage predominant on most transilluminators in the market.

- Compact size with 21 x 26 cm viewing area
- Available in single or dual wavelength models

Cat. No.	Description	List Price	Special Price
U1001	Enduro UV Transilluminator with 302 nm wavelength, 115V	\$1,880.49	<b>\$940.25</b>
U1002	Enduro UV Transilluminator with 302 nm and 365 nm wavelengths, 115V	\$1,985.52	<b>\$893.48</b>

## Labnet Labpette™ FX Fixed Volume Pipettors

Labnet Labpette™ FX Fixed Volume Pipettors are ideal for applications where specific volumes are used repeatedly and sterilization by autoclaving is required.

Cat. No.	Description	List Price	Special Price
P3900-5	5 µL Labpette FX fixed volume pipettor	\$112.97	<b>\$79.08</b>
P3900-10	10 µL Labpette FX fixed volume pipettor	\$114.10	<b>\$79.87</b>
P3900-20	20 µL Labpette FX fixed volume pipettor	\$114.10	<b>\$79.87</b>
P3900-50	50 µL Labpette FX fixed volume pipettor	\$115.23	<b>\$80.66</b>
P3900-250	250 µL Labpette FX fixed volume pipettor	\$115.23	<b>\$80.66</b>
P3900-300	300 µL Labpette FX fixed volume pipettor	\$115.23	<b>\$80.66</b>
P3900-500	500 µL Labpette FX fixed volume pipettor	\$116.38	<b>\$81.47</b>



2LV-7.24

## Labnet BioPette™ Plus Pipettors

BioPette Plus Pipettors feature smooth plunger movement with low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button and is protected by a locking mechanism.

- ▶ Calibration key is included with each pipettor for easy in-house calibration
- ▶ Fully autoclavable and UV-resistant
- ▶ Designed for universal fit with all common brands of pipet tips

Cat. No.	Description	List Price	Special Price
P3942-2	0.1-2.0 µL Biopette Plus pipettor, single-channel	\$273.87	<b>\$191.71</b>
P3942-10	0.5-10 µL Biopette Plus pipettor, single-channel	\$273.87	<b>\$191.71</b>
P3942-20	2-20 µL Biopette Plus pipettor, single-channel	\$257.90	<b>\$180.53</b>
P3942-100	10-100 µL Biopette Plus pipettor, single-channel	\$257.90	<b>\$180.53</b>
P3942-200	20-200 µL Biopette Plus pipettor, single-channel	\$257.90	<b>\$180.53</b>
P3942-1000	100-1,000 µL Biopette Plus pipettor, single-channel	\$257.90	<b>\$180.53</b>
P3942-SK4	Labnet BioPette Plus starter kit	\$933.91	<b>\$653.74</b>



## Labnet Excel™ Electronic Pipettors

The ergonomic, intuitive electronic pipettor for use in most laboratory applications. Its high precision stepping motor controls piston action, resulting in exceptional accuracy and precision.



2LV-9.24

Cat. No.	Description	List Price	Special Price
P3600L-10	0.5-10 µL Excel electronic pipettor, 1-channel	\$394.29	<b>\$276.00</b>
P3600L-20	2-20 µL Excel electronic pipettor, 1-channel	\$394.29	<b>\$276.00</b>
P3600L-200	20-200 µL Excel electronic pipettor, 1-channel	\$394.29	<b>\$276.00</b>
P3600L-1200	100-1,200 µL Excel electronic pipettor, 1-channel	\$394.29	<b>\$276.00</b>
P3608L-10	0.5-10 µL Excel electronic pipettor, 8-channel	\$805.12	<b>\$563.58</b>
P3608L-20	2-20 µL Excel electronic pipettor, 8-channel	\$805.12	<b>\$563.58</b>
P3608L-200	20-200 µL Excel electronic pipettor, 8-channel	\$805.12	<b>\$563.58</b>
P3608L-1200	100-1,200 µL Excel electronic pipettor, 8-channel	\$912.70	<b>\$638.89</b>
P3612L-10	0.5-10 µL Excel electronic pipettor, 12-channel	\$939.79	<b>\$657.85</b>
P3612L-20	2-20 µL Excel electronic pipettor, 12-channel	\$939.79	<b>\$657.85</b>
P3612L-200	20-200 µL Excel electronic pipettor, 12-channel	\$939.79	<b>\$657.85</b>
P3612L-1200	100-1,200 µL Excel electronic pipettor, 12-channel	\$1073.88	<b>\$751.72</b>

**SAVE  
BIG**

Special offers valid  
July 1 – December 31, 2024