



## Standard Series Chromatography Refrigerator

### ABT-12CB

*Temperature Range 0° - 10° C.  
Factory Preset at 4° C*

*1 Swinging Glass Door*

*Exterior Dimensions: 25" W x 23" F-B x 63" H*

*115V, 60Hz, 5 amps  
1/3 HP compressor*

*120 VAC Operation*

*3 Adjustable Shelves*

*Mechanical Thermostat*

*Self-Closing Door*

*Leveling legs*

*CFC Free foam insulation and refrigerant*

*Energy efficient, double pane vacuum insulated glass door*

*Shipping Weight: 275 lbs.*

*Two each 2" access ports with spring-loaded, self-closing gasketed metal cover*

*One each interior GFI protected duplex electrical outlet*



## ABS Standard Chromatography Refrigerator

### ABT-19CB

<p><i>Temperature Range 0° - 10° C. Factory Preset at 4° C</i></p> <p><i>1 Swinging Glass Door</i></p> <p><i>Exterior Dimensions: 27" W x 25" F-B x 75" H</i></p> <p><i>115V, 60Hz, 5 amps 1/3 HP compressor</i></p> <p><i>120 VAC Operation</i></p> <p><i>4 Adjustable Shelves</i></p> <p><i>Mechanical Thermostat</i></p> <p><i>Self-Closing Door</i></p> <p><i>Leveling legs</i></p>	<p><i>CFC Free foam insulation and refrigerant</i></p> <p><i>Energy efficient, double pane vacuum insulated glass door</i></p> <p><i>Shipping Weight: 315 lbs.</i></p> <p><i>Two each 2" access ports with spring-loaded, self-closing gasketed metal cover</i></p> <p><i>One each interior GFI protected duplex electrical outlet</i></p> <p><i>2 Year Parts and Labor Warranty with an additional 3 Year Compressor Parts only warranty</i></p>
---	---