

DOT Mini Vortex Mixer

Description

The Mini vortex mixer is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point. The 3000rpm fixed speed ensures there is no sacrificing any mixing power.

Features

- Extremely compact and easy operation.
- Powerful vortexing for tubes up to 30mm.
- Head is made of durable material to ensure long life.
- Refined steel base is designed to ensure product stability.

Specification

Product Name	Mini Vortex Mixers
Oscillation Mode	Circumference Oscillation
Orbital Diameter	Φ4.5mm
Capacity	Max.80ml
Tube Diameter	≤Φ30
Ambient Temperature	5~40℃
Speed	3000rpm
Power Supply	AC100V-240V / 50Hz-60Hz
Maximum Power	5W
Outside Dimensions (Φ×H)	96×85mm
Weight	0.55kg
IP rating (DIN EN 60529)	IP40