

## Thermo Scientific Owl EC-105 Compact Power Supply

**Ideal for DNA and RNA electrophoresis**



- Equipped with an over current detector for maximum safety
- Easy-to-use power supply operates in constant voltage
- The adjustment knob and hi-lo range switch control the output, which is variable in 10V increments from 20V to 250V
- Two-position switch toggles the bright LED display from indicating output voltage to output current
- The supply is also equipped with an over current detector for maximum safety and has two jacks

**Warranty\*:** 1 year  
**Certifications:** CE marked

Cat. No.	Maximum Voltage	Max. Current	No. of Jacks	Display Modes	L x W x H	Electrical
105ECA-115	250V	300mA @ 150-250V; 500mA below 150V	1	Voltage or current	15 x 20 x 11cm (6 x 8 x 4.5in.)	115V/60Hz
105ECA-LVD	250V	300mA @ 150-250V; 500mA below 150V	1	Voltage or current	15 x 20 x 11cm (6 x 8 x 4.5in.)	230V/50Hz

Cat. No.	Description
112ECA	250V Power Supply Timer

## Thermo Scientific Owl EC-300XL2 Compact Power Supply

**Simple-to-use, multi-purpose, lightweight and great for DNA/RNA and Pprotein electrophoresis**



- Power-off memory retains settings after shut-down
- Stackable design saves benchtop space
- Soft-touch key pad allows quick set-up
- Designed for safety
- Non-skid rubber feet provide stability

**Warranty\*:** 1 year  
**Certifications:** CE marked

Cat. No.	Output Voltage	Max. Current	No. of Jacks	Display Modes	Timer	L x W x H	Electrical
EC300XL	300V	400mA	3	Voltage or current	0-999 min.	26 x 14 x 13.3cm (10.25 x 5.5 x 5.25in.)	120V 50/60Hz
EC300XL2	300V	400mA	3	Voltage or current	0-999 min.	26 x 14 x 13.3cm (10.25 x 5.5 x 5.25in.)	230V 50/60Hz

## Thermo Scientific Owl EC-200XL High Current Power Supply

**Simple-to-use, multi-purpose, lightweight and great for DNA/RNA, protein and blotting electrophoresis**



- Power-off memory retains settings after shut-down
- Stackable design saves benchtop space
- Soft-touch key pad allows quick set-up
- Designed for safety
- Non-skid rubber feet provide stability

**Warranty<sup>Y</sup>:** 1 year  
**Certifications:** CE marked

Cat. No.	Output Voltage	Max. Current	No. of Jacks	Display Modes	Timer	L × W × H	Electrical
EC200XL	200V	2000mA	4	Voltage or current	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	120V 50/60Hz
EC200XL2	200V	2000mA	4	Voltage or current	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	230V 50/60Hz

## Thermo Scientific Owl EC-1000XL Programmable Power Supply

**Great for DNA/RNA, protein and blotting electrophoresis**



- Power-off memory retains settings after shut-down
- Offers constant power mode
- Stackable design saves benchtop space
- Soft-touch key pad allows quick set-up
- Designed for safety
- Non-skid rubber feet
- Easy to operate

**Warranty<sup>Y</sup>:** 1 year  
**Certifications:** CE marked

Cat. No.	Output Voltage	Max. Current	No. of Jacks	Display Modes	Timer	L × W × H	Electrical
EC1000XL	1000V	500mA	4	Voltage, current, or power	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	120V 50/60Hz
EC1000XL2	1000V	500mA	4	Voltage, current, or power	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	230V 50/60Hz

## Thermo Scientific Owl EC-3000XL Programmable Power Supply

**Simple-to-use, multi-purpose and lightweight**



- Power-off memory retains settings after shut-down
- Offers constant power mode
- Stackable design saves benchtop space
- Soft-touch key pad allows for quick set-up
- Designed for safety
- Non-skid rubber feet

**Warranty<sup>Y</sup>:** 1 year  
**Certifications:** CE marked

Cat. No.	Output Voltage	Max. Current	No. of Jacks	Display Modes	Timer	L × W × H	Electrical
<b>EC3000XL</b>	3000V	400mA	4	Voltage, current, or power	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	120V 50/60Hz
<b>EC3000XL2</b>	3000V	400mA	4	Voltage, current, or power	0-99 hr. 59 min.	31 × 27 × 13.3cm (12 × 10.75 × 5.25in.)	230V 50/60Hz

### Thermo Scientific Owl Power Supply Leads

- Power supply leads include a retractable sheath end for convenience and safety
- Power supply leads are manufactured with strain relief cord to prevent breakage

Cat. No.	Description
<b>PSL-5</b>	Power supply leads, 3 ft. (includes one black and one red lead)