


LABNET INTERNATIONAL, INC.

FASTPETTE^{PLUS}

Pipette Controller

Labnet
Labnet International, Inc.




- **Ergonomic design for exceptional comfort**
- **Digital LCD display for continuous visualization of pipetting speed and battery status**
- **Rechargeable battery with convenient charging stand**
- **High aspiration and dispensing speed, 25ml in under 3.5 seconds**



FASTPETTE PLUS™

The New **FASTPETTE Plus** is a high performance pipette controller designed for cordless operation with glass or plastic pipettes from 0.5 to 100ml. With a sleek, ergonomic handle design and contoured push buttons, the controller fits comfortably in the hand even after long periods of use. The included charging stand can be used for storing the pipettor, at the same time, charging the battery.

The large LCD display indicates the battery status along with the pipetting mode and speed. Aspiration/dispensing speeds can be set to high or low to adjust for the volume of the pipette or the viscosity of the sample. The high setting provides quick filling of large volume pipettes and can aspirate as much as 25ml in 3.2 seconds. The low speed setting is ideal for precise measurement of smaller volumes.



Ordering Information

Item No.	Description:
P2001	FASTPETTE Plus Pipette Controller with charging stand, Universal Plug*
P2023	PTFE autoclavable filter, 0.2µm, pk/5
P2024	PTFE autoclavable filter, 0.45µm, pk/5
P2028	Wall hanger (included)
P2001-CS	Charging stand (included)

* Universal power supply is suitable for US, EU, UK and AU plug configurations.