

# Accessories

## SPECIFICATIONS TABLE/ORDER NUMBERS


DESCRIPTION	ORDER NUMBER	DETAILS
<b>Additional shelving – gravity convection</b>		
Additional wire mesh shelf, 60 L models	50127761	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 12.28 x 15.20 in
Additional wire mesh shelf, 100 L models	50127762	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 16.61 x 15.20 in
Additional wire mesh shelf, 180 L models	50127763	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 16.61 x 22.07 in
Additional wire mesh shelf, 400 L models	50135245	Wire mesh shelf for General Protocol ovens; incl. 2 shelf supports; 20.79 x 21.06 in
Additional wire mesh shelf, 750 L models	50135246	Wire mesh shelf for General Protocol ovens; incl. 2 shelf supports; 38.74 x 21.06 in
Reinforced shelf, 180 L models	50128880	Shelf with maximum capacity of 32 lb for 180 L models, gravity convection of Advanced Protocol and Advanced Protocol Security ovens; max total capacity is 114 lb
<b>Additional shelving – mechanical convection</b>		
Additional wire mesh shelf, 60 L models	50127764	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 13.31 x 13.22 in
Additional wire mesh shelf, 100 L models	50127765	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 13.64 x 13.22 in
Additional wire mesh shelf, 180 L models	50127766	Wire mesh shelf for General Protocol; Advanced Protocol; Advanced Protocol Security ovens; incl. 2 shelf supports; 17.64 x 20.12 in
Additional wire mesh shelf, 400 L models	50135245	Wire mesh shelf for Advanced Protocol ovens; incl. 2 shelf supports; 20.79 x 21.06 in
Additional wire mesh shelf, 750 L models	50135246	Wire mesh shelf for Advanced Protocol ovens; incl. 2 shelf supports; 38.74 x 21.06 in
Stainless steel perforated shelf, 60 L models	50127773	Stainless steel perforated shelf for General Protocol oven; incl. 2 shelf supports; 12.95 x 14.80 in
Stainless steel perforated shelf, 100 L models	50127774	Stainless steel perforated shelf for General Protocol oven; incl. 2 shelf supports; 17.28 x 14.80 in
Stainless steel perforated shelf, 180 L models	50127777	Stainless steel perforated shelf for General Protocol oven; incl. 2 shelf supports; 17.28 x 21.69 in
Reinforced shelf, 180 L models	50128881	Shelf with maximum capacity of 32 lb for 180 L models, mechanical convection of Advanced Protocol and Advanced Protocol Security ovens; max total capacity is 114 lb
<b>Sample sensor</b>		
Sample sensor, for select models	50127767	Sample sensor for connection to Advanced Protocol Security table top ovens (60 L, 100 L, 180 L); large capacity Advanced Protocol large capacity ovens (400 L, 750 L): measures exact sample temperature; sample temperature is shown on display as plugged in; cable length: 86.61 in
<b>Silicone free viton door sealing</b>		
Silicone free viton door sealing, 60 L models	50130657	Silicone free viton door sealing for all 60 L ovens
Silicone free viton door sealing, 100 L models	50130658	Silicone free viton door sealing for all 100 L ovens
Silicone free viton door sealing, 180 L models	50130659	Silicone free viton door sealing for all 180 L ovens
Silicone free viton door sealing, 400 L models	50135869	Silicone free viton door sealing for all 400 L ovens
Silicone free viton door sealing, 750 L models	50135870	Silicone free viton door sealing for all 750 L ovens
<b>Fresh air filter</b>		
Fresh air particle filter, all ovens models	50127566	Fresh air particle filter for connection to the fresh air vent (intake); for all ovens models
<b>Stacking kits</b>		
Stacking kit, all 60 L models	50126665	Stacking kit to stack two 60 L models
Stacking kit, all 100 L models	50126666	Stacking kit to stack two 100 L models, or 60 L on 100 L
Stacking kit, all 180 L models	50126667	Stacking kit to stack two 100 L models, or 60 L /100 L on 180 L
<b>Support stands</b>		
Support stand, all 60 L models	50127741	Support stand with casters for 60 L models; height including casters: 7.36 in
Support stand, all 100 L models	50127742	Support stand with casters for 100 L models; height including casters: 7.36 in
Support stand, all 180 L models	50127743	Support stand with casters for 180 L models; height including casters: 7.36 in



Stacking kit



Support stand with casters



Wire mesh shelves

a. Stainless steel perforated shelves  
b. Reinforced shelves



Fresh air particle filter