


TECHNOLOGY

EBA 21

Small centrifuge, without rotor		
Power supply ^{*)}	220–240 V 1 ~	110–127 V 1 ~
Frequency	50–60 Hz	
Consumption	500 VA	550 VA
Emission	EN 55011 group 1, class B EN 61000-3-2, EN 61000-3-3	FCC class B
Immunity	EN 61000-6-2	-
 Max. capacity		
	6x50 ml	
Max. RPM (speed)		
	18,000 min ⁻¹	
Max. RCF		
	25,718	
Running time	1 sec – 99 min : 59 sec, ∞ continuous run, short cycle mode (impulse key)	
Dimensions (H x W x D)		
	247 x 275 x 330 mm	
Weight		
	approx. 11 kg	
Cat. No.	1004	1004-01

^{*)} Other voltages on request.



Hettich centrifuges comply with all relevant EU standards in effect and conform to the European level of quality and safety for medical devices. Evidence is provided by national and international test marks such as IEC 61010 or the CE conformity. The ISO 9001:2008 and ISO 13485:2003 certificates accredited to the company bear witness to the extreme care and responsibility Hettich puts into the manufacturing of centrifuges and their accessories.

Andreas Hettich GmbH & Co. KG

Föhrenstr. 12
D - 78532 Tuttlingen
Germany

www.hettichlab.com
info@hettichlab.com
service@hettichlab.com

Phone +49 (0)7461 / 705 -0
Fax +49 (0)7461 / 705 -122

National Sales: -200
International Sales: -201
National Service: -202
International Service: -203