




# ROTANTA 460/460 R

TECHNOLOGY	ROTANTA 460		ROTANTA 460 R <sup>1&gt;)</sup>		
<b>Benchtop centrifuge, without rotor</b>	<b>classic</b>		<b>cooled</b>		
Power supply	200 – 240 V 1 ~	100 – 127 V 1 ~	200 – 240 V 1 ~	100–127 V 1~	100 V 1 ~
Frequency	50 – 60 Hz		50 Hz	60 Hz	50 Hz
Consumption	1,000 VA	1,100 VA	1,800 VA	2,000 VA	
Emission	EN / IEC 61326-1, class B	FCC class B	EN / IEC 61326-1, class B	FCC class B	
Immunity		-		-	
 Max. capacity	4x750 ml				
 Max. RPM (speed)	15,000 min <sup>-1</sup>				
 max. RCF	24,400				
Running time	1 sec - 99 h : 59 min : 59 sec, ∞ continuous run, short cycle mode				
Dimensions (HxWxD)	456 x 554 x 706 mm		456 x 770 x 706 mm		
Weight	approx. 93 kg	approx. 103 kg	approx. 133 kg	approx. 143 kg	
<b>Refrigeration</b>					
Temperature control, infinitely variable	-		from -20 to +40 °C		
<b>Cat. No.</b>	<b>5650</b>	<b>5650-01</b>	<b>5660</b>	<b>5660-01</b>	

	ROTANTA 460 RC	ROTANTA 460 RF
	<b>Underbench centrifuge, cooled, without rotor</b>	<b>Floorstanding centrifuge, cooled, without rotor</b>
Power supply	200 – 240 V 1 ~	
Frequency	50 Hz	
Consumption	1,800 VA	
Emission / Immunity	EN / IEC 61326-1, class B	
Dimensions (HxWxD)	683 x 554 x 697 mm	961 x 554 x 697 mm
Weight	approx. 140 kg	approx. 164 kg
Temperature control, infinitely variable	from -20 to +40 °C	
	<b>The other values are identical with those of the ROTANTA 460 R</b>	
<b>Cat. No.</b>	<b>5670</b>	<b>5675</b>

- \*) Other voltages on request.  
 >) Additional model variants on request.



Hettich centrifuges comply with all relevant EU standards in effect and conform to the European level of quality and safety for medical devices. Evidence is provided by national and international test marks such as IEC 61010 or the CE conformity. The ISO 9001:2008, ISO 13485:2003 and 14001:2004 certificates accredited to the company bear witness to the extreme care and responsibility Hettich puts into the manufacturing of centrifuges and their accessories.

## Andreas Hettich GmbH & Co. KG

Föhrenstr. 12  
 D - 78532 Tuttlingen  
 Germany

www.hettichlab.com  
 info@hettichlab.com  
 service@hettichlab.com

Phone +49 (0)7461 / 705 -0  
 Fax +49 (0)7461 / 705 -122

**National Sales: -200**  
**International Sales: -201**  
**National Service: -202**  
**International Service: -203**

