



Product Specifications

20X SSC BUFFER CONCENTRATE

<u>Catalog #</u>	<u>Size</u>
IB72010	1L

Physical Specifications

CAS#	N/A
Solution Components	3.0M Sodium Chloride 0.3M Sodium Citrate
Conductivity (1X Solution)	17,750 ± 2150µmhos
pH (20X Solution @ 25°C)	7.00 ± .15
Refractive Index (20X Solution)	1.370 ± 0.010

Molecular Biology Specifications

DNase assay	None Detected
RNase assay	None Detected

Recommended Use

20X SSC is a solution formulated for use in nucleic acid preparation and blotting applications. It is used in concentrations ranging from 0.2X to 20X, depending on the applications.

Storage

Store at room temperature

Warning

Irritant. Avoid contact with eyes, skin and mucous membranes. Avoid breathing vapors. Wash thoroughly with water after handling. See Material Safety Data Sheet for additional information.