

Smoothie™ 1 ml–25ml Pipette Fillers

The newest addition to our range of pipette fillers, the Jencons Smoothie™ is exceptionally smooth in action, simple to use and flexible in application, eliminating the need to use a range of different sized pipette fillers for a variety of volumes dispensed.

You can choose from a range of five fabulous colors each impregnated with its own subtle fragrance: Vanilla, Banana, Blueberry, Strawberry or Lemon & Lime*.

- Comfortable and simple to use
- Smooth and controllable thumb screw action
- Precise and accurate pipetting
- Compatible with blow-out, TD (to deliver) or gravity dispense pipettes
- Fits pipettes from 1–25ml
- Removable and autoclavable silicone pipette holder
- Acrylic stand available accommodating up to five Smoothies
- 3 year warranty

The Smoothie's unique range selector allows you to select from two dispense resolutions for maximum accuracy. One range setting delivers fast filling with optimal accuracy when used with 10ml to 25ml pipettes, the other setting is for use with pipettes from 1ml to 10ml where accuracy of fill can sometimes be difficult to control.

Cat. No.	Description	Color
475-360	Vanilla Smoothie™ pipette filler	White and blue
475-361	Strawberry Smoothie™ pipette filler	Red and green
475-362	Banana Smoothie™ pipette filler	Yellow and gray
475-363	Lemon & Lime Smoothie™ pipette filler	Green and yellow
475-364	Blueberry Smoothie™ pipette filler	Blue and purple
475-002	Stand for five Smoothie™ pipette fillers	Clear acrylic
266-004	Replacement silicone pipette holder	

* Over a period of time the fragrance will gradually dissipate.



Safetypette™ Pipette Filler

A simple yet controllable pipette filler, the Safetypette™ is designed to give superior pipetting control to left and right-handed users. The lightweight Safetypette™ sits comfortably in the hand, enabling smooth control in both aspirating and dispensing, using 1ml to 100ml pipettes. Simply squeeze the silicone bulb and the thumb lever controls both the aspirate and dispense modes. Incorporates central button to blow out residual liquid in blow-out pipettes. Supplied with three additional colored pipette cones.

- Comfortable and simple to use
- Precise pipetting control
- Robust and lightweight
- Large capacity bulb
- Compatible with blow-out pipettes
- Simple thumb button operation
- Uses replaceable standard membrane filter
- Autoclavable silicone pipette holder
- Supplied with integral 0.45µm filter and colored pipette cones for safety color coding in white, mint, mauve and pink
- Easily maintained and cleaned
- Three year warranty



Button to expel residual liquid when using blow-out pipettes.



Comfortable to use



Supplied with integral but replaceable 0.45µm filter



Cat. No. Description

Cat. No.	Description	Pack
475-102	Safetypette™ pipette filler with integral 0.45µm filter and colored pipette cones set	–
266-004	Replacement silicone pipette holder	1
266-003	Replacement filter set 0.45µm	5
266-089	Replacement filter set 0.2µm	5
266-090	Replacement sterile filter 0.2µm	1

Safetypette™ Pro Pipette Filler

The Safetypette™ Pro is a unique design of pipette filler with an ultra-squeezable bulb to provide smooth, manual pipetting control in both aspirating and dispensing, using 1ml to 100ml pipettes.

Simply squeeze the large silicone bulb and the thumb lever controls both the aspirate and dispense modes, with a button to blow-out residual contents of pipette if required.

- Comfortable and simple to use
- Precise pipetting control
- Robust and lightweight
- Large capacity bulb
- Compatible with blow-out pipettes
- Simple thumb button operation
- Uses integral 0.45µm, replaceable membrane filter to ensure liquid is not accidentally drawn into the unit
- Autoclavable silicone pipette holder
- Three year warranty



Button to expel residual liquid when using blow-out pipettes



Comfortable to use



Supplied with integral but replaceable 0.45µm filter



Cat. No.	Description	Pack
475-155	Safetypette™ Pro pipette filler with integral 0.45µm filter	1
266-004	Replacement silicone pipette holder	1
266-003	Replacement filter set 0.45µm	5
266-089	Replacement filter set 0.2µm	5
266-090	Replacement sterile filter 0.2µm	1