

# MEGAMAG GENIE®

## Large Volume Magnetic Stirrers



**ANALOG**



**DIGITAL**



**LOW SPEED**

**LOW RANGE**  
**LR**  
**LOW SPEED**  
**VERSION NOW**  
**AVAILABLE**

**High Speed Analog and Digital or Low Speed Digital - Stir up to 38 Liter Containers!**

**Another revolutionary product from the makers of the Vortex-Genie®!**

- ⦿ Choice of THREE different models - Analog or Digital for high speed applications, and Low Range Digital for low speed applications.
- ⦿ Powerful "rare earth" magnets stay coupled even under the most demanding stirring conditions.
- ⦿ Rugged ABS/Polycarbonate housing and the famous Vortex-Genie motor.
- ⦿ Digital models offer a smooth ramp-up to desired speed avoiding potential spills.
- ⦿ Digital models incorporate a timer from 1 to 99 minutes.
- ⦿ May be used in cold rooms or incubators.



### Ordering Information

Catalog No.	Description
SI-2236	MegaMag Genie, 120V
SI-2246	MegaMag Genie, 230V, No Plug
SI-2256	MegaMag Genie, 230V, Euro Plug
SI-2266	MegaMag Genie, 230V, British Plug
SI-2276	MegaMag Genie, 230V, Swiss Plug
SI-2296	MegaMag Genie, 240V, Australian Plug
SI-2286	MegaMag Genie, 100V
SI-3236	Digital MegaMag Genie, 120V
SI-3246	Digital MegaMag Genie, 230V, No Plug
SI-3256	Digital MegaMag Genie, 230V, Euro Plug
SI-3266	Digital MegaMag Genie, 230V, British Plug
SI-3276	Digital MegaMag Genie, 230V, Swiss Plug
SI-3296	Digital MegaMag Genie, 240V, Australian Plug
SI-3286	Digital MegaMag Genie, 100V
SI-3236L	Digital MegaMag Genie LR, 120V
SI-3246L	Digital MegaMag Genie LR, 230V, No Plug
SI-3256L	Digital MegaMag Genie LR, 230V, Euro Plug
SI-3266L	Digital MegaMag Genie LR, 230V, British Plug
SI-3276L	Digital MegaMag Genie LR, 230V, Swiss Plug
SI-3296L	Digital MegaMag Genie LR, 240V, Australian Plug
SI-3286L	Digital MegaMag Genie LR, 100V

### Specifications

Top Surface Dimensions (D x W)	305 x 305mm (12 x 12in)
Overall Dimensions (D x W x H)	310 x 305 x 130mm (12.25 x 12 x 5.1in)
Weight	3 kgs (6.6 lbs)
Stirring Speed	250-1300 RPM 50-300 RPM (LR model) (50-285 for 50Hz)
Time	1-99 min or continuous (digital models only)
Maximum Stirring Positions	1
Maximum Stirring Capacity	38 Liters
Stir Bar Maximum Length	7.6cm (3 in)
Material	ABS/Polycarbonate



Home of the Multi-Task Genie™ Family  
70 Orville Drive, Bohemia, NY 11716 U.S.A.  
Phone: (631) 567-4700 • Fax: (631) 567-5896  
Toll Free: 888-850-6208  
Customer Service: [custsvic@scientificindustries.com](mailto:custsvic@scientificindustries.com)  
[scientificindustries.com](http://scientificindustries.com)