

Description

The 3D nutating shaker product is an economical mini 3D shaker for mixing samples with the combined motions of an orbital shaker and a rocker. This provides a gentle yet thorough mixing of the sample.

Features

- Small footprint and an elegant design.
- 3D motion generates thorough and gentle mixing.
- Dimpled non-skid mat holds tubes firmly in place.

Specification

Model	NS-01A	Previous Model	BC-4201A
Speed	20rpm	Platform Dimensions	260×260mm
Tilt Angle	20°	Dimensions	260×260×190mm
Max. Load	0.8Kg	Power	20W
Ambient Temperature	5°C~40°C	Power Supply	AC100V/60Hz(US) AC220V/50Hz(Eur)
Relative Humidity	≤80%	Weight	1.8Kg

