



PV-1 personal vortex mixer

Extremely compact personal vortex mixer with a low profile and small footprint for gentle mixing to vigorous re-suspension of cell or chemical pellets in 0.5 ml, 1.5 ml and 15 ml tubes.

- Adjustable speed control: 750 to 3000 rpm
- 'Continuous' or 'touch' operation
- For tubes up to 20 mm diameter

Pressure sensitive cup accommodates tubes up to 20 mm diameter

Reliable and extremely quiet motor produces regulated and reproducible agitation throughout the speed range

Extremely easy to operate – select either 'continuous' or 'touch' operation and turn the knob to adjust the speed from 750 rpm to 3000 rpm

Low voltage cord easily fits through incubator door gaskets

In 'touch' mode, agitation starts in response to pressure on the pressure sensitive cup

Exceptionally compact with a low profile and small footprint – fits into almost any location

Rubber suction pads hold tight to the work surface and prevent the unit from 'walking' – they also absorb vibration and prevent its transmission to the workbench

