

Options (Field-installed must be installed by a qualified professional)

Description	For use with ▶	UxF300	UxF400	UxF500	UxF600	UxF700
LN₂ Backup System Maintains temperature down to -80°C with liquid nitrogen Models for UxF300 are stand alone.	Factory Installed	N/A	LN4567	LN4567	LN4567	LN4567
	Field Installed	FLN300A* FLN300D* FLN300V* FLN300W*	FLN4567	FLN4567	FLN4567	FLN4567
CO₂ Backup System Maintains temperature down to -67°C with CO ₂ Models for UxF300 are stand alone.	Factory Installed	N/A	CO4567	CO4567	CO4567	CO4567
	Field Installed	FCO300A* FCO300D* FCO300V* FCO300W*	FCO4567	FCO4567	FCO4567	FCO4567
Chart Recorder 6" (15.2cm), seven day inkless recorder, to -150°C	Factory Installed	N/A	CR4567	CR4567	CR4567	CR4567
	Field Installed	6283-6A* 6283-6D* 6283-6V* 6283-6W*	FCR400R	FCR567R	FCR567R	FCR567R
Access Key Option Card access control Includes five key cards.	Factory Installed	RAC34567	RAC34567	RAC34567	RAC34567	RAC34567
	Field Installed	FRAC34567	FRAC34567	FRAC34567	FRAC34567	FRAC34567














* A = 115V/60 Hz (NEMA 5-20 plug), D = 208-230V/60 Hz (NEMA 6-15 plug), V = 230V/50 Hz (EU plug), W = 230V/50 Hz (British plug)

Accessories

Description	For use with ▶	UxF300	UxF400	UxF500	UxF600	UxF700
Access Key Pack (U.S.) Includes five cards; supports ISO15693 Protocol		ACU34567	ACU34567	ACU34567	ACU34567	ACU34567
Access Key Pack (EU) Includes five cards; supports ISO14443 Protocol		ACE34567	ACE34567	ACE34567	ACE34567	ACE34567
Shelf Kit Includes one shelf and clips		SK300	SK400	SK500	SK600	SK700
Chart Paper inkless (pack of 50)		N/A	6185	6185	6185	6185
Replacement Air Filter		AF34567	AF34567	AF34567	AF34567	AF34567
Replacement Back-Up Battery		400159	400159	400159	400159	400159
Alarm Delay Module Eliminates nuisance alarms due to intermittent or transitory conditions. Adjustable delay requires alarm condition to exist for user defined period before signal is released to monitoring system.		6903	6903	6903	6903	6903
Cryo Gloves® Medium		4425	4425	4425	4425	4425
Cryo Gloves Large		4426	4426	4426	4426	4426

Standard Plug Option Configurations

Choose one of our specialty plugs to configure your freezer to meet local requirements

Model No.	Plug Description	Male End
Standard Plugs		
US125V20A	NEMA 5-20	
US230V15A	NEMA 6-15	
EU230V16A	Continental EU	
Specialty Plugs		
IT230V16A	Italy	
UK230V13A	Great Britain	
IN230V16A	India	
IS230V16A	Israel	
AU230V16A	Australia	
CH230V16A	China	
DK230V16A	Denmark	
AR230V16A	Argentina	
SW230V16A	Switzerland	
BR230V16A	Brazil	



Cryo Gloves



Inkless Chart Paper





Access Key Option








Shelf Kit

Thermo Scientific Revco UxF Series

Thermo Scientific Revco UxF Series Ultra-Low Temperature Upright Freezer Racks

2" Box Racks	Model No. Description	Dimensions H x W x D in. (cm)	Storage	UxF300	UxF400	UxF500	UxF600	UxF700
	398329 Sliding Drawer Rack for 2" Boxes	9.5 x 5.5 x 26.3 (24 x 14 x 66.8)	Boxes Per Rack	20	20	20	20	20
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Boxes Per Freezer	300	400	500	600	700
	1950518 Side Access Rack for 2" Boxes	9.4 x 5.4 x 26.5 (23.9 x 13.7 x 69.9)	Boxes Per Rack	20	20	20	20	20
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Boxes Per Freezer	300	400	500	600	700

3" Box Racks	Model No. Description	Dimensions H x W x D in. (cm)	Storage	UxF300	UxF400	UxF500	UxF600	UxF700
	398328 Sliding Drawer Rack for 3" Boxes	9.5 x 5.5 x 26.3 (24 x 14 x 66.8)	Boxes Per Rack	15	15	15	15	15
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Boxes Per Freezer	225	300	375	450	525
	1950519 Side Access Rack for 3" Boxes	9.4 x 5.4 x 26.5 (23.9 x 13.7 x 67.3)	Boxes Per Rack	15	15	15	15	15
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Boxes Per Freezer	225	300	375	450	525

Microplate Racks	Model No. Description	Dimensions H x W x D in. (cm)	Storage	UxF300	UxF400	UxF500	UxF600	UxF700
	1950639 Sliding Drawer Rack for Standard or Deepwell	9.5 x 5.5 x 26.3 (24 x 14 x 66.8)	Plates Per Rack	28	28	28	28	28
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Plates Per Freezer	420	560	700	840	980
	1950522 Side Access Rack for Standard Plates	9.4 x 5.5 x 25.5 (23.9 x 14 x 64.8)	Plates Per Rack	77	77	77	77	77
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Plates Per Freezer	1155	1540	1925	2310	2695
	398322 Side Access Rack with Locking Rod for Standard or Deepwell	9.4 x 5.5 x 25.7 (23.9 x 14 x 65.3)	Plates Per Rack	98	98	98	98	98
			Racks Per Shelf	3	4	5	6	7
			Racks Per Freezer	15	20	25	30	35
			Plates Per Freezer	294	392	490	588	686

Racking Shelf Kits Includes Sliding Drawer Racks and 2" Boxes with 100 Count Cell Dividers						
		UxF300	UxF400	UxF500	UxF600	UxF700
	Model No. ▶	RSK300SD5	RSK400SD5	RSK500SD5	RSK600SD5	RSK700SD5
	Racks Included	3	4	5	6	7.0
	Boxes Included	60	80	100	120	140