

S2030-LS Orbit™ LS Low Speed Shaker



The Orbit™ LS low speed shaker is a high quality, general purpose shaker with economical analog control.

- Economical analog control with easy-turn knobs
- 19 mm orbit for blotting, gel staining and general mixing
- Streamlined design
- Safe for cold room and incubator use

Styled after the other models in the line, the Orbit™ LS features a streamlined, compact design with easy-access control panel. Electronic controls are utilized for setting speed and time. The knobs used for selecting parameters turn easily and are clearly marked with corresponding values.

A large orbit and low speed make the Orbit™ LS useful for staining and destaining gels, washing blots and general mixing. The shaker is supplied with a 30 x 30 cm stainless steel flat platform and non-slip rubber mat to keep plates, dishes and other sample containers in place.

Safe for use in temperature controlled environments, the small profile of this shaker makes it ideal for use in incubators such as the Labnet 311 and 611 series.

PRODUCT SPECIFICATIONS	
Speed range	3 - 60 rpm
Timer	20 - 120 min (10 min inc.) or continuous
Motion	Circular, 19 mm
Maximum load	11 lb/5 kg
Ambient operating range	+4° to 70°C
Dimensions (W x D x H)	10 x 12.6 x 6.3 in/25.5 x 32 x16 cm
Weight	14.5 lb/6.6 kg
Electrical	230V~ or 120V~, 50/60 Hz

