

Transilluminators



Standard Series, Single Wavelength UV

- 15-watt units provide high UV intensity
- Fixed-intensity models include TS-, TC-, TR- and TL-series
- Variable-intensity models with UV-blocking cover include TV-series
- Models accommodate a variety of sample sizes
- Housing dimensions (W x L x H): 18.5" x 19.25" x 4", 34.3 x 48.9 x 10.2cm
- Available in 120V/60Hz, 230V/50Hz, 240V/50Hz and 100V/50-60Hz

Our high power (15-watt) **Standard Series** units provide unsurpassed UV intensities and UV irradiance to produce brilliant fluorescent responses and highest clarity. Ideal for revealing even trace amounts of DNA! Models available in various UV wavelengths, fixed and variable intensity, filter sizes, timer options and UV-blocking covers.

Fixed-Intensity TS-, TC-, TR- and TL-Series

Model	Wavelength	Filter Size	Filter Assembly	Tubes	Features
TS-254R	254nm	6" x 6", 15 x 15cm	2F606B	(6) 15W SW, BLE-1T155	S
TS-312R	312nm	6" x 6", 15 x 15cm	3F606R	(6) 15W MW, BLE-1T158	S, D
TS-365R	365nm	6" x 6", 15 x 15cm	2F326B	(6) 15W LW, BLE-1T151	S, D
TC-254R	254nm	8" x 8", 20 x 20cm	2F808B	(6) 15W SW, BLE-1T155	T
TC-312R	312nm	8" x 8", 20 x 20cm	3F808R	(6) 15W MW, BLE-1T158	T, D
TC-365R	365nm	8" x 8", 20 x 20cm	2F883B	(6) 15W LW, BLE-1T151	T, D
TR-254R	254nm	6" x 14", 15 x 35cm	2F614B	(6) 15W SW, BLE-1T155	S
TR-312R	312nm	6" x 14", 15 x 35cm	2F614B	(6) 15W MW, BLE-1T158	S, D
TR-365R	365nm	9" x 13", 23 x 33cm	2F365B	(6) 15W LW, BLE-1T151	S, D
TL-254R	254nm	8" x 16", 20 x 40cm	2F816B	(6) 15W SW, BLE-1T155	T
TL-312R	312nm	8" x 16", 20 x 40cm	2F816B	(6) 15W MW, BLE-1T158	T, D
TL-365R	365nm	8" x 16", 20 x 40cm	2F168B	(6) 15W LW, BLE-1T151	T, D

Variable-Intensity TV-Series

Model	Wavelength	Filter Size	Filter Assembly	Tubes	Features
TVC-312R	312nm	8" x 8", 20 x 20cm	3F808R	(6) 15W MW, BLE-1T158	S, D, C, V
TVR-312R	312nm	6" x 14", 15 x 35cm	2F614B	(6) 15W MW, BLE-1T158	S, D, C, V
TVL-312R	312nm	8" x 16", 20 x 40cm	2F816B	(6) 15W MW, BLE-1T158	S, D, C, V

