

MATERIAL SAFETY DATA SHEET



IBI SCIENTIFIC
9861 KAPP COURT
PEOSTA, IA 52068

EMERGENCY TELEPHONE NUMBER:
3E COMPANY (800) 451-8346
(001) 760-602-8703

TELEPHONE NUMBER FOR INFORMATION:
(800) 253-4942 OR (563) 690-0484
E-MAIL: INFO@IBISCI.COM

DATE PREPARED: 1-12-99

DATE REVISED: 1-14-10

SECTION 1. ----- CHEMICAL IDENTIFICATION-----

CATALOG#: IB70054 & IB70055
NAME: AGAROSE, ULTRA SIEVE
CAS#: 9012-36-6

SECTION 2. ----- HAZARDS IDENTIFICATION-----

LABEL PRECAUTIONARY STATEMENT

MAY BE IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.

IN CASE OF CONTACT WITH EYES, RINSE IMMEDIATELY WITH PLENTY OF
WATER AND SEEK MEDICAL ADVICE.

WEAR SUITABLE PROTECTIVE CLOTHING.

CAUTION: THE TOXOLOGICAL PROPERTIES OF THIS CHEMICAL HAVE NOT BEEN
THOROUGHLY INVESTIGATED.

HMS: NFPA:

H: 1

H: 1

F: 0

F: 0

R: 1

R: 0

PE: B

SH: N/A

SECTION 3. ----- COMPOSITION/INFORMATION ON INGREDIENTS -----

THIS MATERIAL IS NOT CONSIDERED HAZARDOUS UNDER 29CFR 1910.1200. AS WITH
ALL CHEMICALS EXERCISE CAUTION WHILE HANDLING.

SECTION 4. ----- FIRST-AID MEASURES-----

IN CASE OF CONTACT, IMMEDIATELY FLUSH EYES WITH COPIOUS AMOUNTS OF
WATER FOR AT LEAST 15 MINUTES.

IN CASE OF CONTACT, IMMEDIATELY WASH SKIN WITH SOAP AND COPIOUS
AMOUNTS OF WATER.

IF INHALED, REMOVE TO FRESH AIR. IF NOT BREATHING GIVE ARTIFICIAL
RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN.

IF SWALLOWED, WASH OUT MOUTH WITH WATER PROVIDED PERSON IS CONSCIOUS.
CALL A PHYSICIAN.

WASH CONTAMINATED CLOTHING BEFORE REUSE.

SECTION 5. ----- FIRE FIGHTING MEASURES -----

EXTINGUISHING MEDIA

WATER SPRAY.

CARBON DIOXIDE, DRY CHEMICAL POWDER OR APPROPRIATE FOAM.

SPECIAL FIREFIGHTING PROCEDURES

WEAR SELF-CONTAINED BREATHING APPARATUS AND PROTECTIVE CLOTHING TO PREVENT CONTACT WITH SKIN AND EYES.

UNUSUAL FIRE AND EXPLOSIONS HAZARDS

EMITS TOXIC FUMES UNDER FIRE CONDITIONS.

SECTION 6. - - - - - ACCIDENTAL RELEASE MEASURES - - - - -

WEAR RESPIRATOR, CHEMICAL SAFETY GOGGLES, RUBBER BOOTS AND HEAVY RUBBER GLOVES.

SWEEP UP, PLACE IN A BAG AND HOLD FOR WASTE DISPOSAL.

AVOID RAISING DUST.

VENTILATE AREA AND WASH SPILL SITE AFTER MATERIAL PICKUP IS COMPLETE.

SECTION 7. - - - - - HANDLING AND STORAGE - - - - -

DO NOT INHALE.

AVOID CONTACT WITH EYES, SKIN AND CLOTHING.

WASH THOROUGHLY AFTER HANDLING.

KEEP TIGHTLY CLOSED.

STORE IN A COOL DRY PLACE.

STORE AT ROOM TEMPERATURE.

SECTION 8. - - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION - - - - -

WEAR APPROPRIATE NIOSH/MSHA-APPROVED RESPIRATOR, CHEMICAL-RESISTANT GLOVES, SAFETY GOGGLES, OTHER PROTECTIVE CLOTHING.

SAFETY SHOWER AND EYE BATH.

CHEMICAL FUME HOOD RECOMMENDED.

SECTION 9. - - - - - PHYSICAL AND CHEMICAL PROPERTIES - - - - -

APPEARANCE AND ODOR

SOLID.

SECTION 10. - - - - - STABILITY AND REACTIVITY - - - - -

STABILITY

STABLE.

INCOMPATIBILITIES

STRONG OXIDIZERS

HAZARDOUS COMBUSTION OR DECOMPOSITION PRODUCTS

TOXIC FUMES OF:

CARBON MONOXIDE, CARBON DIOXIDE. VETILATE AREA AND WASH SPIL SITE AFTER MATERIAL PICKUP IS COMPLETE.

HAZARDOUS POLYMERIZATION

WILL NOT OCCUR.

SECTION 11. - - - - - TOXICOLOGICAL INFORMATION - - - - -

ACUTE EFFECTS

MAY BE HARMFUL IF SWALLOWED.

DUST MAY BE IRRITATING TO THE EYES, MUCOUS MEMBRANES AND UPPER RESPIRATORY TRACT.

MAY CAUSE SKIN IRRITATION.

TO THE BEST OF OUR KNOWLEDGE, THE CHEMICAL, PHYSICAL, AND

TOXICOLOGICAL PROPERTIES HAVE NOT BEEN THOROUGHLY INVESTIGATED.

SECTION 12. ECOLOGICAL INFORMATION

DATA NOT YET AVAILABLE.

SECTION 13. DISPOSAL CONSIDERATIONS

OBSERVE ALL FEDERAL, STATE AND LOCAL ENVIRONMENTAL REGULATIONS.

SECTION 14. TRANSPORT INFORMATION

NOT REGULATED.

SECTION 15. REGULATORY INFORMATION

EUROPEAN INFORMATION

R: N/A

S: N/A

SECTION 16. OTHER INFORMATION.....

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT BUT DOES NOT PURPORT TO BE ALL INCLUSIVE AND SHALL BE USED ONLY AS A GUIDE. IBI SCIENTIFIC SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR FROM CONTACT WITH THE ABOVE PRODUCT. SEE REVERSE SIDE OF INVOICE OR PACKING SLIP FOR ADDITIONAL TERMS AND CONDITIONS OF SALE.