

ProBlot™ Rocker 35 - Adjustable Tilt Rocker

- ❑ *Single and double platform models*
- ❑ *Durable construction with stainless steel platforms*
- ❑ *Variable tilt angle for a wide variety of applications*
- ❑ *Ideal for staining gels and membranes*
- ❑ *Safe for cold room or incubator use*

Specifications:

Speed range	4-160 tilts/min
Timer	0-120 min. or continuous
Motion	rocking
Tilt angle	adjustable +/-20°
Max. load	7kg
Dimensions (WxDxH)	14x11.4x6.5in/36x29x17cm
Platform dimensions	14x11in/36x28cm
Double platform clearance	3.5in./9cm
Ambient operating range	4° to 65°C

Ordering Information:

S2035	<i>ProBlot 35 Rocker with platform, 115V</i>
S2035-230V	<i>ProBlot 35 Rocker with platform, 230V</i>
S2035-D	<i>ProBlot 35D Rocker with double platform, 115V</i>
S2035-D-230V	<i>ProBlot 35D Rocker with double platform, 230V</i>
S2025-ST	<i>Staining trays, pack of 4</i>



S2035-D

